

Product datasheet for TL502437

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

Yap1 Mouse shRNA Plasmid (Locus ID 22601)

Product data:

Product Type: shRNA Plasmids

Product Name: Yap1 Mouse shRNA Plasmid (Locus ID 22601)

Locus ID: 2260°

Synonyms: Al325207; Y; Yap; Yap65; Yk; Yki; yor; Yorkie

Vector: pGFP-C-shLenti (TR30023)

E. coli Selection: Chloramphenicol (34 ug/ml)

Mammalian Cell

Selection:

Puromycin

Format: Lentiviral plasmids

Components: Yap1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 22601).

5µg purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.

RefSeq: BC014733, BC039125, BC094313, NM 001171147, NM 009534, NM 009534.1, NM 009534.2,

NM 009534.3, NM 001171147.1

UniProt ID: P46938

Summary: This gene encodes a protein which binds to the SH3 domain of the Yes proto-oncogene

product, a tyrosine kinase. This protein contains a WW domain, consisting of four conserved aromatic amino acids including two tryptophan residues. This conserved WW domain is found in various structural, regulatory and signaling molecules in various species, and may play a role in protein-protein interaction. Following cellular damage, phosphorylation of this

encoded protein may suppress apoptosis. This protein may be involved in malignant

transformation in cancer. Alternative splicing results in multiple transcript variants. [provided

by RefSeq, Jan 2010]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>.

If you need a special design or shRNA sequence, please utilize our custom shRNA service.





Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).