

Product datasheet for **TL502609**

Plod3 Mouse shRNA Plasmid (Locus ID 26433)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Plod3 Mouse shRNA Plasmid (Locus ID 26433)
Locus ID:	26433
Synonyms:	AI414586; LH3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Plod3 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 26433). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC043047 , BC054734 , NM_011962 , NM_011962.1 , NM_011962.2 , NM_011962.3
UniProt ID:	Q9R0E1
Summary:	Multifunctional enzyme that catalyzes a series of post-translational modifications on Lys residues in procollagen (PubMed:16447251). Plays a redundant role in catalyzing the formation of hydroxylysine residues in -Xaa-Lys-Gly- sequences in collagens (PubMed:16447251). Plays a redundant role in catalyzing the transfer of galactose onto hydroxylysine groups, giving rise to galactosyl 5-hydroxylysine (By similarity). Has an essential role by catalyzing the subsequent transfer of glucose moieties, giving rise to 1,2-glucosylgalactosyl-5-hydroxylysine residues (PubMed:16447251, PubMed:16467571, PubMed:21220425). Catalyzes hydroxylation and glycosylation of Lys residues in the MBL1 collagen-like domain, giving rise to hydroxylysine and 1,2-glucosylgalactosyl-5-hydroxylysine residues (PubMed:25419660). Catalyzes hydroxylation and glycosylation of Lys residues in the ADIPOQ collagen-like domain, giving rise to hydroxylysine and 1,2-glucosylgalactosyl-5-hydroxylysine residues (PubMed:23209641). Essential for normal biosynthesis and secretion of type IV collagens (PubMed:15377789, PubMed:16467571, PubMed:17873278). Essential for normal formation of basement membranes (PubMed:15377789, PubMed:16467571). [UniProtKB/Swiss-Prot Function]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).