

## Product datasheet for **TL502651V**

### PPP2R5C Mouse shRNA Lentiviral Particle (Locus ID 26931)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	PPP2R5C Mouse shRNA Lentiviral Particle (Locus ID 26931)
Locus ID:	26931
Synonyms:	2610043M05Rik; 2700063L20Rik; AI060890; AW545884; C85228; D12Bwg0916e; mKIAA0044
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	PPP2R5C - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">BC003979</a> , <a href="#">NM_001081457</a> , <a href="#">NM_001081458</a> , <a href="#">NM_001135001</a> , <a href="#">NM_012023</a> , <a href="#">NM_001081458.1</a> , <a href="#">NM_001081458.2</a> , <a href="#">NM_012023.1</a> , <a href="#">NM_012023.2</a> , <a href="#">NM_012023.3</a> , <a href="#">NM_001081457.1</a> , <a href="#">NM_001081457.2</a> , <a href="#">NM_001135001.1</a> , <a href="#">BC147179</a> , <a href="#">BC147180</a> , <a href="#">NM_001364204</a> , <a href="#">NM_001364205</a> , <a href="#">NM_001364206</a> , <a href="#">NM_001364207</a> , <a href="#">NM_001364208</a> , <a href="#">NM_001364209</a> , <a href="#">NM_001364210</a> , <a href="#">NM_001364211</a> , <a href="#">NM_001364212</a> , <a href="#">NM_001364213</a> , <a href="#">NM_001364214</a> , <a href="#">NM_001081458.3</a> , <a href="#">NM_001081457.3</a> , <a href="#">NM_001135001.2</a> , <a href="#">NM_012023.4</a>
UniProt ID:	<a href="#">Q60996</a>
Summary:	The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment. The PP2A-PPP2R5C holoenzyme may activate TP53 and play a role in DNA damage-induced inhibition of cell proliferation. PP2A-PPP2R5C may also regulate the ERK signaling pathway through ERK dephosphorylation (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).