

Product datasheet for **TL502790**

Zbtb18 Mouse shRNA Plasmid (Locus ID 30928)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Zbtb18 Mouse shRNA Plasmid (Locus ID 30928)
Locus ID:	30928
Synonyms:	RP58; zfp-238; Zfp238; Znf238
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Zbtb18 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 30928). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC054529 , BC054742 , NM_001012330 , NM_013915 , NM_001355098 , NM_001355099 , NM_001355100 , NM_001355101 , NM_001355102 , NM_001355103 , NM_013915.1 , NM_013915.2 , NM_013915.3 , NM_001012330.1 , BC050909 , BC089333
UniProt ID:	Q9WUK6
Summary:	Transcriptional repressor that plays a role in various developmental processes such as myogenesis and brain development. Specifically binds the consensus DNA sequence 5'-[AC]ACATCTG[GT][AC]-3' which contains the E box core, and acts by recruiting chromatin remodeling multiprotein complexes. Plays a key role in myogenesis by directly repressing the expression of ID2 and ID3, 2 inhibitors of skeletal myogenesis. Also involved in controlling cell division of progenitor cells and regulating the survival of postmitotic cortical neurons. May also play a role in the organization of chromosomes in the nucleus.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).