

## Product datasheet for **TL502863**

### Atg13 Mouse shRNA Plasmid (Locus ID 51897)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Atg13 Mouse shRNA Plasmid (Locus ID 51897)
Locus ID:	51897
Synonyms:	1110053A20Rik; D2Ertd391e; Harbi1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Atg13 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 51897). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC016669</a> , <a href="#">BC023673</a> , <a href="#">BC023712</a> , <a href="#">BC033419</a> , <a href="#">NM_145528</a> , <a href="#">NM_001355416</a> , <a href="#">NM_001355417</a> , <a href="#">NM_001355418</a> , <a href="#">NM_001355419</a> , <a href="#">NM_001355420</a> , <a href="#">NM_001355421</a> , <a href="#">NM_145528.1</a> , <a href="#">NM_145528.2</a> , <a href="#">NM_145528.3</a> , <a href="#">BC024457</a>
UniProt ID:	<a href="#">Q91YI1</a>
Summary:	Autophagy factor required for autophagosome formation and mitophagy. Target of the TOR kinase signaling pathway that regulates autophagy through the control of the phosphorylation status of ATG13 and ULK1, and the regulation of the ATG13-ULK1-RB1CC1 complex. Through its regulation of ULK1 activity, plays a role in the regulation of the kinase activity of mTORC1 and cell proliferation.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).