

## Product datasheet for **TL502902**

### Suz12 Mouse shRNA Plasmid (Locus ID 52615)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Suz12 Mouse shRNA Plasmid (Locus ID 52615)
Locus ID:	52615
Synonyms:	2610028O16Rik; AI195385; AU016842; AW536442; D11Ert530; D11Ert530e; mKIAA0160
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Suz12 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 52615). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC064461</a> , <a href="#">NM_001163018</a> , <a href="#">NM_199196</a> , <a href="#">NM_199196.1</a> , <a href="#">NM_199196.2</a> , <a href="#">NM_001163018.1</a> , <a href="#">BC003922</a> , <a href="#">BC031375</a> , <a href="#">BC039915</a> , <a href="#">BC051099</a> , <a href="#">BC084591</a>
UniProt ID:	<a href="#">Q80U70</a>
Summary:	This gene encodes a core component of the polycomb repressive complex 2 (PRC2) that also includes, at least, embryonic ectoderm development protein (EED) and enhancer of zeste homolog 1 or 2 (EZH1 or EZH2). Through the methyltransferase activity of EZH1 or EZH2, the PRC2 complex methylates Lys9 and Lys27 of histone 3 and Lys26 of histone 1, leading to recruitment of the PRC1 complex, histone 2A ubiquitylation and transcriptional repression of the target genes. This gene product is essential for the activity and integrity of the PRC2 complex, and is required for X chromosome inactivation, stem cell maintenance and differentiation. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).