

Product datasheet for **TL503055**

Atp8b1 Mouse shRNA Plasmid (Locus ID 54670)

Product data:

Product Type:	shRNA Plasmids
Locus ID:	54670
Synonyms:	AI451886; FIC1; Ic
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Atp8b1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector (Gene ID = 54670). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<u>BC117946</u> , <u>NM_001001488</u> , <u>NM_001001488.1</u> , <u>NM_001001488.2</u> , <u>NM_001001488.3</u>
UniProt ID:	<u>Q148W0</u>
Summary:	Catalytic component of a P4-ATPase flippase complex which catalyzes the hydrolysis of ATP coupled to the transport of aminophospholipids from the outer to the inner leaflet of various membranes and ensures the maintenance of asymmetric distribution of phospholipids. Phospholipid translocation seems also to be implicated in vesicle formation and in uptake of lipid signaling molecules. May play a role in asymmetric distribution of phospholipids in the canicular membrane. Plays a role in bile salt homeostasis. In cooperation with ABCB4 may be involved in establishing integrity of the canicular membrane thus protecting hepatocytes from bile salts. Involved in the microvillus formation in polarized epithelial cells; the function seems to be independent from its flippase activity. Required for the preservation of cochlear hair cells in the inner ear. Required for the preservation of cochlear hair cells in the inner ear. According PubMed:20852622 is proposed to act as cardiolipin transporter during inflammatory injury; the function is questioned by PubMed:21475228.[UniProtKB/Swiss-Prot Function]



shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).