

Product datasheet for TL503108V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Abcb10 Mouse shRNA Lentiviral Particle (Locus ID 56199)

Product data:

Product Type: shRNA Lentiviral Particles

Product Name: Abcb10 Mouse shRNA Lentiviral Particle (Locus ID 56199)

Locus ID: 56199

Synonyms: ABC-m; Abc-me; Abcb12

Vector: pGFP-C-shLenti (TR30023)

Format: Lentiviral particles

Components: Abcb10 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble

control), 0.5 ml each, >10^7 TU/ml.

RefSeq: <u>BC046818</u>, <u>BC053020</u>, <u>BC054793</u>, <u>NM 019552</u>, <u>NM 019552.1</u>, <u>NM 019552.2</u>, <u>BC053020.1</u>

UniProt ID: Q9JI39

Summary: This gene encodes a member of the ATP-binding cassette superfamily of transporters. ATP-

binding cassette proteins transport various molecules across extra- and intra-cellular

membranes. The encoded protein is localized to the mitochondrial inner membrane where it interacts with and stabilizes mitoferrin-1, and is important for heme biosynthesis. Additional

evidence suggests the encoded protein is involved in oxidative stress protection and

erythropoisesis. [provided by RefSeq, May 2013]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our custom shRNA service.





Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).