

Product datasheet for **TL503139**

Arl6 Mouse shRNA Plasmid (Locus ID 56297)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Arl6 Mouse shRNA Plasmid (Locus ID 56297)
Locus ID:	56297
Synonyms:	1110018H24Rik; 2210411E14Rik; BBS3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Arl6 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 56297). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC018497 , NM_001347244 , NM_019665 , NM_019665.1 , NM_019665.2 , NM_019665.3
UniProt ID:	O88848
Summary:	Involved in membrane protein trafficking at the base of the ciliary organelle (By similarity). Mediates recruitment onto plasma membrane of the BBSome complex which would constitute a coat complex required for sorting of specific membrane proteins to the primary cilia (By similarity). Together with BBS1, is necessary for correct trafficking of PKD1 to primary cilia (PubMed:24939912). Together with the BBSome complex and LTZL1, controls SMO ciliary trafficking and contributes to the sonic hedgehog (SHH) pathway regulation (By similarity). May regulate cilia assembly and disassembly and subsequent ciliary signaling events such as the Wnt signaling cascade (By similarity). Isoform 2 may be required for proper retinal function and organization (PubMed:20333246).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).