

## Product datasheet for **TL503192**

### Mtch2 Mouse shRNA Plasmid (Locus ID 56428)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Mtch2 Mouse shRNA Plasmid (Locus ID 56428)
Locus ID:	56428
Synonyms:	2310034D24Rik; 4930539J07Rik; HSPC0; HSPC032
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Mtch2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 56428). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC038899</a> , <a href="#">BC038916</a> , <a href="#">NM_019758</a> , <a href="#">NM_019758.1</a> , <a href="#">NM_019758.2</a> , <a href="#">NM_019758.3</a> , <a href="#">BC002152</a> , <a href="#">BC028929</a> , <a href="#">BC065145</a>
UniProt ID:	<a href="#">Q791V5</a>
Summary:	This gene encodes a member of the SLC25 family of nuclear-encoded transporters that are localized in the inner mitochondrial membrane. Members of this superfamily are involved in many metabolic pathways and cell functions. Genome-wide association studies in human have identified single-nucleotide polymorphisms in several loci associated with obesity. This gene is one such locus, which is highly expressed in white adipose tissue and adipocytes, and thought to play a regulatory role in adipocyte differentiation and biology. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. A recent study showed this gene to be an authentic stop codon readthrough target that can produce two isoforms from the same mRNA by use of alternative in-frame translation termination codons. [provided by RefSeq, Dec 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).