

Product datasheet for **TL503354V**

Cnot9 Mouse shRNA Lentiviral Particle (Locus ID 58184)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cnot9 Mouse shRNA Lentiviral Particle (Locus ID 58184)
Locus ID:	58184
Synonyms:	2610007F23Rik; AI593551; F110; Rqcd1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Rqcd1 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC050898 , BC051948 , NM_021383 , NM_021383.1 , NM_021383.2 , NM_021383.3 , NM_021383.4 , NM_021383.5 , BC034663
UniProt ID:	Q9JKY0
Summary:	Component of the CCR4-NOT complex which is one of the major cellular mRNA deadenylases and is linked to various cellular processes including bulk mRNA degradation, miRNA-mediated repression, translational repression during translational initiation and general transcription regulation. Additional complex functions may be a consequence of its influence on mRNA expression. Involved in down-regulation of MYB- and JUN-dependent transcription. May play a role in cell differentiation. Required for retinoic acid-induced differentiation of F9 teratocarcinoma cells. Does not bind DNA by itself. Enhances ligand-dependent transcriptional activity of nuclear hormone receptors. May play a role in cell differentiation. [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).