

Product datasheet for **TL503388**

Wrap73 Mouse shRNA Plasmid (Locus ID 59002)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Wrap73 Mouse shRNA Plasmid (Locus ID 59002)
Locus ID:	59002
Synonyms:	2610044M17Rik; 5330425N03Rik; Dd57; Wdr8; Wrap73
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Wrap73 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 59002). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC016120 , NM_021499 , NM_021499.1 , NM_021499.2 , BC064111
UniProt ID:	Q9JM98
Summary:	The SSX2IP:WRAP73 complex is proposed to act as regulator of spindle anchoring at the mitotic centrosome. Required for the centrosomal localization of SSX2IP and normal mitotic bipolar spindle morphology. Required for the targeting of centriole satellite proteins to centrosomes such as of PCM1, SSX2IP, CEP290 and PIBF1/CEP90. Required for ciliogenesis and involved in the removal of the CEP97:CCP110 complex from the mother centriole. Involved in ciliary vesicle formation at the mother centriole and required for the docking of vesicles to the basal body during ciliogenesis; may promote docking of RAB8A- and ARL13B-containing vesicles (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).