

Product datasheet for **TL503489**

Nif3l1 Mouse shRNA Plasmid (Locus ID 65102)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Nif3l1 Mouse shRNA Plasmid (Locus ID 65102)
Locus ID:	65102
Synonyms:	1110030G24Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Nif3l1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 65102). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC065163 , BC089475 , NM_022988 , NM_022988.1 , NM_022988.2 , NM_022988.3 , BC129866
UniProt ID:	Q9EQ80
Summary:	This gene is a member of the NGG1-interacting factor 3-like superfamily of transcriptional regulators and is ubiquitously expressed throughout embryonic development. The encoded protein interacts with a component of the constitutive photomorphogenesis 9 signalosome, and functions as a transcriptional corepressor of genes involved in neuronal differentiation. This gene is highly conserved from bacteria to human. [provided by RefSeq, Sep 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).