

## Product datasheet for **TL503738**

### Nacc1 Mouse shRNA Plasmid (Locus ID 66830)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Nacc1 Mouse shRNA Plasmid (Locus ID 66830)
Locus ID:	66830
Synonyms:	2010001H03Rik; 4930511N13Rik; Btbd14b; Nac1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Nacc1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 66830). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC048157</a> , <a href="#">BC052706</a> , <a href="#">BC054800</a> , <a href="#">NM_025788</a> , <a href="#">NM_025788.1</a> , <a href="#">NM_025788.2</a> , <a href="#">NM_025788.3</a>
UniProt ID:	<a href="#">Q7TSZ8</a>
Summary:	Functions as a transcriptional repressor. Seems to function as a transcriptional corepressor in neuronal cells through recruitment of HDAC3 and HDAC4. Contributes to tumor progression, and tumor cell proliferation and survival. This may be mediated at least in part through repressing transcriptional activity of GADD45GIP1. Required for recruiting the proteasome from the nucleus to the cytoplasm and dendritic spines. Involved in the acute behavioral and neurological responses to cocaine and amphetamines.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).