

Product datasheet for **TL503800**

Fancl Mouse shRNA Plasmid (Locus ID 67030)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Fancl Mouse shRNA Plasmid (Locus ID 67030)
Locus ID:	67030
Synonyms:	2010322C19Rik; AW554273; B230118H11Rik; gcd; P; Phf; Phf9; Pog
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Fancl - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 67030). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC032876 , BC066181 , NM_001277273 , NM_025923 , NR_102382 , NM_025923.1 , NM_025923.2 , NM_025923.3 , NM_001277273.1 , BC032876.1
UniProt ID:	Q9CR14
Summary:	This gene encodes the complementation group L subunit of the multimeric Fanconi anemia (FA) nuclear complex composed of proteins encoded by over ten Fanconi anemia complementation (FANC) group genes. The FA complex is necessary for protection against DNA damage. This gene product, an E3 ubiquitin ligase, catalyzes and is required for the monoubiquitination of the protein encoded by the Fanconi anemia, complementation group D2 gene, a critical step in the FA pathway (PMID: 12973351, 21229326). In mouse, mutations of this E3 ubiquitin ligase gene can lead to infertility in adult males and females, and a deletion of this gene can cause embryonic lethality in some genetic backgrounds. A pseudogene of this gene has been identified on chromosome 1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).