

## Product datasheet for **TL503947**

### Tbc1d23 Mouse shRNA Plasmid (Locus ID 67581)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Tbc1d23 Mouse shRNA Plasmid (Locus ID 67581)
Locus ID:	67581
Synonyms:	4930451A13Rik; AU015720; AU043671; AU043778; D030022P07Rik; W51689
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Tbc1d23 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 67581). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC031706</a> , <a href="#">NM_026254</a> , <a href="#">NM_001358433</a> , <a href="#">NM_001358434</a> , <a href="#">NM_026254.1</a> , <a href="#">NM_026254.2</a> , <a href="#">BC019762</a>
UniProt ID:	<a href="#">Q8K0F1</a>
Summary:	Putative Rab GTPase-activating protein which plays a role in vesicular trafficking. Involved in endosome-to-Golgi trafficking. Acts as a bridging protein by binding simultaneously to golgins, including GOLGA1 and GOLGA4, located at the trans-Golgi, and to the WASH complex, located on endosome-derived vesicles (PubMed:29084197). Together with WDR11 complex facilitates the golgin-mediated capture of vesicles generated using AP-1 (By similarity). Plays a role in brain development, including in cortical neuron positioning. May also be important for neurite outgrowth, possibly through its involvement in membrane trafficking and cargo delivery, 2 processes which are essential for axonal and dendritic growth (PubMed:28823707). May act as a general inhibitor of innate immunity signaling, strongly inhibiting multiple TLR and dectin/CLEC7A-signaling pathways. Does not alter initial activation events, but instead affects maintenance of inflammatory gene expression several hours after bacterial lipopolysaccharide (LPS) challenge (PubMed:22312129). [UniProtKB/Swiss-Prot Function]



[View online »](#)

**shRNA Design:**

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).