

Product datasheet for **TL504363**

Zfp219 Mouse shRNA Plasmid (Locus ID 69890)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Zfp219 Mouse shRNA Plasmid (Locus ID 69890)
Locus ID:	69890
Synonyms:	2010302A17Rik; Znf219
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Zfp219 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 69890). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC071271 , BC096486 , NM_001253694 , NM_001253695 , NM_001253696 , NM_027248 , NM_027248.1 , NM_027248.2 , NM_001253696.1 , NM_001253694.1 , NM_001253695.1 , BC010447 , BC034874 , BC052050
UniProt ID:	Q6IQX8
Summary:	Transcriptional regulator (PubMed:20940257). Recognizes and binds 2 copies of the core DNA sequence motif 5'-GGGGG-3' (PubMed:20940257). Binds to the HMGN1 promoter and may repress HMGN1 expression (By similarity). Regulates SNCA expression in primary cortical neurons (By similarity). Binds to the COL2A1 promoter and activates COL2A1 expression, as part of a complex with SOX9 (PubMed:20940257). Plays a role in chondrocyte differentiation (PubMed:20940257).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).