

## Product datasheet for **TL504556**

### **Bbs7 Mouse shRNA Plasmid (Locus ID 71492)**

#### **Product data:**

|                                  |   |
|----------------------------------|---|
| <b>Product Type:</b>             | shRNA Plasmids  |
| <b>Product Name:</b>             | Bbs7 Mouse shRNA Plasmid (Locus ID 71492)   |
| <b>Locus ID:</b>                 | 71492   |
| <b>Synonyms:</b>                 | 8430406N16Rik   |
| <b>Vector:</b>                   | pGFP-C-shLenti (TR30023)  |
| <b>E. coli Selection:</b>        | Chloramphenicol (34 ug/ml)  |
| <b>Mammalian Cell Selection:</b> | Puromycin   |
| <b>Format:</b>                   | Lentiviral plasmids   |
| <b>Components:</b>               | Bbs7 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 71492).<br>5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.   |
| <b>RefSeq:</b>                   | <a href="#">BC031505</a> , <a href="#">NM_027810</a> , <a href="#">NM_027810.1</a> , <a href="#">NM_027810.2</a> , <a href="#">NM_027810.3</a>  |
| <b>UniProt ID:</b>               | <a href="#">Q8K2G4</a>  |
| <b>Summary:</b>                  | The BBSome complex is thought to function as a coat complex required for sorting of specific membrane proteins to the primary cilia. The BBSome complex is required for ciliogenesis but is dispensable for centriolar satellite function. This ciliogenic function is mediated in part by the Rab8 GDP/GTP exchange factor, which localizes to the basal body and contacts the BBSome. Rab8(GTP) enters the primary cilium and promotes extension of the ciliary membrane. Firstly the BBSome associates with the ciliary membrane and binds to RAB3IP/Rabin8, the guanosyl exchange factor (GEF) for Rab8 and then the Rab8-GTP localizes to the cilium and promotes docking and fusion of carrier vesicles to the base of the ciliary membrane. The BBSome complex, together with the LTZL1, controls SMO ciliary trafficking and contributes to the sonic hedgehog (SHH) pathway regulation. Required for BBSome complex ciliary localization but not for the proper complex assembly (By similarity).<br>[UniProtKB/Swiss-Prot Function] |
| <b>shRNA Design:</b>             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .  |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).