

Product datasheet for **TL504887**

Msl1 Mouse shRNA Plasmid (Locus ID 74026)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Msl1 Mouse shRNA Plasmid (Locus ID 74026)
Locus ID:	74026
Synonyms:	2810017F12Rik; 4121402D02Rik; 4930463F05Rik; AA682082; Msl-1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Msl1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 74026). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC043039 , BC055715 , BC058629 , NM_028722 , NM_028722.1 , NM_028722.2 , BM946681 , NM_001361931 , NM_001361932 , NM_001361933 , NM_001361934 , NM_028722.3
UniProt ID:	Q6PDM1
Summary:	Component of histone acetyltransferase complex responsible for the majority of histone H4 acetylation at 'Lys-17' which is implicated in the formation of higher-order chromatin. structure (By similarity). Greatly enhances MSL2 E3 ubiquitin ligase activity, promoting monoubiquitination of histone H2B at 'Lys-35' (H2BK34Ub). This modification in turn stimulates histone H3 methylation at 'Lys-5' (H3K4me) and 'Lys-80' (H3K79me) and leads to gene activation, including that of HOXA9 and MEIS1. In the MSL complex, acts as a scaffold to tether MSL3 and KAT8 together for enzymatic activity regulation.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).