

Product datasheet for **TL504918**

Opa1 Mouse shRNA Plasmid (Locus ID 74143)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Opa1 Mouse shRNA Plasmid (Locus ID 74143)
Locus ID:	74143
Synonyms:	1200011N24Rik; AI225888; AI847218; lilr3; mKIAA0567
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Opa1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 74143). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001199177 , NM_133752 , NM_133752.1 , NM_133752.2 , NM_133752.3 , NM_001199177.1 , BC138665 , BC025160 , BC052704 , BC145959
UniProt ID:	P58281
Summary:	Dynamin-related GTPase that is essential for normal mitochondrial morphology by regulating the equilibrium between mitochondrial fusion and mitochondrial fission (PubMed:11847212, PubMed:24616225, PubMed:26785494, PubMed:28746876). Coexpression of isoform 1 with shorter alternative products is required for optimal activity in promoting mitochondrial fusion (By similarity). Binds lipid membranes enriched in negatively charged phospholipids, such as cardiolipin, and promotes membrane tubulation. The intrinsic GTPase activity is low, and is strongly increased by interaction with lipid membranes (By similarity). Plays a role in remodeling cristae and the release of cytochrome c during apoptosis (PubMed:16839884, PubMed:16839885). Proteolytic processing in response to intrinsic apoptotic signals may lead to disassembly of OPA1 oligomers and release of the caspase activator cytochrome C (CYCS) into the mitochondrial intermembrane space (PubMed:16839884, PubMed:16839885). Plays a role in mitochondrial genome maintenance (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).