

Product datasheet for **TL505312V**

Dist Mouse shRNA Lentiviral Particle (Locus ID 78920)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Dist Mouse shRNA Lentiviral Particle (Locus ID 78920)
Locus ID:	78920
Synonyms:	1600017E01Rik; 4632413C10Rik; 4930529O08Rik; DLTS
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Dist - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC006702 , NM_030225 , NM_030225.1 , NM_030225.2 , NM_030225.3 , NM_030225.4 , BC024066
UniProt ID:	Q9D2G2
Summary:	Dihydrolipoamide succinyltransferase (E2) component of the 2-oxoglutarate dehydrogenase complex (By similarity). The 2-oxoglutarate dehydrogenase complex catalyzes the overall conversion of 2-oxoglutarate to succinyl-CoA and CO(2) (By similarity). The 2-oxoglutarate dehydrogenase complex is mainly active in the mitochondrion. A fraction of the 2-oxoglutarate dehydrogenase complex also localizes in the nucleus and is required for lysine succinylation of histones: associates with KAT2A on chromatin and provides succinyl-CoA to histone succinyltransferase KAT2A (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL505312A virus into HEK293 cells. TL505312A virus was prepared using lenti-shRNA TL505312A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL505312B virus into HEK293 cells. TL505312B virus was prepared using lenti-shRNA TL505312B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL505312C] virus into HEK293 cells. [TL505312C] virus was prepared using lenti-shRNA [TL505312C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL505312D] virus into HEK293 cells. [TL505312D] virus was prepared using lenti-shRNA [TL505312D] and [TR30037] packaging kit.