

Product datasheet for **TL505340**

Pofut2 Mouse shRNA Plasmid (Locus ID 80294)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Pofut2 Mouse shRNA Plasmid (Locus ID 80294)
Locus ID:	80294
Synonyms:	2310011G23Rik; AI256847; BC003494; C21orf80; FUT13
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Pofut2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 80294). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_030262 , NM_030262.1 , NM_030262.3 , BC147078 , BC003494 , BC018194 , BC147077
UniProt ID:	Q8VHI3
Summary:	Catalyzes the reaction that attaches fucose through an O-glycosidic linkage to a conserved serine or threonine residue in the consensus sequence C1-X(2,3)-S/T-C2-X(2)-G of thrombospondin type I repeats (TSRs) where C1 and C2 are the first and second cysteines of the repeat, respectively. O-fucosylates members of several protein families including the ADAMTS superfamily and the thrombosporin (TSP) and spondin families. Required for the proper secretion of ADAMTS family members such as ADAMSL1 and ADAMST13 (By similarity). O-fucosylation of TSRs is also required for restricting epithelial to mesenchymal transition (EMT), maintaining the correct patterning of mesoderm and localization of the definite endoderm.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).