

## Product datasheet for **TL505485**

### Sirt1 Mouse shRNA Plasmid (Locus ID 93759)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Sirt1 Mouse shRNA Plasmid (Locus ID 93759)
Locus ID:	93759
Synonyms:	AA673258; S; Si; Sir2; Sir2a; Sir2alpha; SIR2L1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Sirt1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 93759). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001159589</a> , <a href="#">NM_001159590</a> , <a href="#">NM_019812</a> , <a href="#">NM_019812.1</a> , <a href="#">NM_019812.2</a> , <a href="#">NM_019812.3</a> , <a href="#">NM_001159589.1</a> , <a href="#">NM_001159589.2</a> , <a href="#">NM_001159590.1</a> , <a href="#">BC152314</a> , <a href="#">BC006584</a> , <a href="#">BC165995</a>
UniProt ID:	<a href="#">Q923E4</a>
Summary:	This gene encodes a member of the sirtuin family of proteins, characterized by their deacetylase activity and proposed role in longevity. The encoded protein regulates gene expression in a wide range of cell and tissue types through its NAD <sup>+</sup> -dependent deacetylation of histones, transcription factors and transcriptional coactivators. Brain-specific overexpression of this gene has been shown to result in increased median lifespan. Viability of homozygous knockout mice for this gene varies with strain background. Homozygous knockout mice of strains that do not exhibit embryonic lethality are sterile and have a reduced lifespan. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).