

Product datasheet for **TL505803**

Adgrb1 Mouse shRNA Plasmid (Locus ID 107831)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Adgrb1 Mouse shRNA Plasmid (Locus ID 107831)
Locus ID:	107831
Synonyms:	B830018M07Rik; Bai1; mKIAA4089; R75078
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Adgrb1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 107831). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_174991 , NM_001359759 , NM_174991.1 , NM_174991.2 , NM_174991.3 , BC036963 , BC059857 , BC115835
UniProt ID:	Q3UHD1
Summary:	Phosphatidylserine receptor which enhances the engulfment of apoptotic cells (PubMed:17960134). Also mediates the binding and engulfment of Gram-negative bacteria (PubMed:21245295, PubMed:26838550, PubMed:26838550). Stimulates production of reactive oxygen species by macrophages in response to Gram-negative bacteria, resulting in enhanced microbicidal macrophage activity (By similarity). In the gastric mucosa, required for recognition and engulfment of apoptotic gastric epithelial cells (By similarity). Promotes myoblast fusion (PubMed:23615608). Activates the Rho pathway in a G-protein-dependent manner (By similarity). Inhibits MDM2-mediated ubiquitination and degradation of DLG4/PSD95, promoting DLG4 stability and regulating synaptic plasticity (PubMed:25751059). Required for the formation of dendritic spines by ensuring the correct localization of PARD3 and TIAM1 (By similarity). Potent inhibitor of angiogenesis in brain and may play a significant role as a mediator of the p53/TP53 signal in suppression of glioblastoma (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>