

## Product datasheet for **TL505834V**

### Taf9 Mouse shRNA Lentiviral Particle (Locus ID 108143)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Taf9 Mouse shRNA Lentiviral Particle (Locus ID 108143)
Locus ID:	108143
Synonyms:	2310012M09Rik; AA673500; AU045423; Taf2g; TAFII31; TAFII32
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Taf9 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">BC019453</a> , <a href="#">BC042723</a> , <a href="#">BC043028</a> , <a href="#">BC108348</a> , <a href="#">NM_001015889</a> , <a href="#">NM_027139</a> , <a href="#">NM_027592.1</a> , <a href="#">NM_027592.2</a> , <a href="#">NM_027139.1</a> , <a href="#">NM_027139.2</a> , <a href="#">NM_027139.3</a> , <a href="#">NM_027139.4</a> , <a href="#">NM_027139.5</a> , <a href="#">NM_001015889.1</a> , <a href="#">NM_001015889.2</a> , <a href="#">BC043028.1</a> , <a href="#">BC096683</a>
UniProt ID:	<a href="#">Q8VI33</a>
Summary:	Essential for cell viability. TAF9 and TAF9B are involved in transcriptional activation as well as repression of distinct but overlapping sets of genes. May have a role in gene regulation associated with apoptosis. TAFs are components of the transcription factor IID (TFIID) complex, the TBP-free TAFII complex (TFTC), the PCAF histone acetylase complex and the STAGA transcription coactivator-HAT complex. TFIID or TFTC are essential for the regulation of RNA polymerase II-mediated transcription (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).