

## Product datasheet for **TL505910**

### Zfp91 Mouse shRNA Plasmid (Locus ID 109910)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Zfp91 Mouse shRNA Plasmid (Locus ID 109910)
Locus ID:	109910
Synonyms:	Pzf, AL024263, AW545902, 9130014I08Rik, A530054C17Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Zfp91 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 109910). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001039718</a> , <a href="#">NM_053009</a> , <a href="#">NM_053009.1</a> , <a href="#">NM_053009.2</a> , <a href="#">NM_053009.3</a> , <a href="#">BC057323</a> , <a href="#">BC064766</a> , <a href="#">BC083000</a> , <a href="#">BC156276</a> , <a href="#">BC157098</a> , <a href="#">BC168394</a>
UniProt ID:	<a href="#">Q62511</a>
Summary:	The protein encoded by this gene is a member of the zinc finger family of proteins. The gene product contains C2H2-type domains, which are the classical zinc finger domains found in numerous nucleic acid-binding proteins. The homologous human protein has been shown to function as a regulator of the non-canonical NF-kappaB pathway in lymphotoxin-beta receptor signaling. A read-through transcript variant composed of Zfp91 and the downstream Cntf gene sequence has been identified, but it is thought to be non-coding. Read-through transcription of Zfp91 and Cntf has been observed in both human and mouse. A Zfp91-related pseudogene has also been identified on chromosome 17. [provided by RefSeq, Oct 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).