

Product datasheet for **TL505926**

Triobp Mouse shRNA Plasmid (Locus ID 110253)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Triobp Mouse shRNA Plasmid (Locus ID 110253)
Locus ID:	110253
Synonyms:	AI428493; EST478828; mKIAA1662; Tara
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Triobp - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 110253). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001024716 , NM_001039155 , NM_001039156 , NM_138579 , NM_001024716.1 , NM_001039156.1 , NM_138579.1 , NM_138579.2 , NM_138579.3 , NM_138579.4 , NM_001039155.1 , BC003984 , BC042760 , BC053329 , BC094366
UniProt ID:	Q99KW3
Summary:	This gene encodes a protein that interacts with trio, which is involved with neural tissue development and controlling actin cytoskeleton organization, cell motility, and cell growth. The encoded protein also associates with F-actin and stabilizes F-actin structures. Domains contained in this encoded protein are an N-terminal pleckstrin homology domain and a C-terminal coiled-coil region. Mutations in the human gene have been associated with a form of autosomal recessive nonsyndromic deafness. Multiple alternatively spliced transcript variants have been described [provided by RefSeq, Sep 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).