

Product datasheet for **TL506295**

Mtmr9 Mouse shRNA Plasmid (Locus ID 210376)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Mtmr9 Mouse shRNA Plasmid (Locus ID 210376)
Locus ID:	210376
Synonyms:	9430075G12Rik; AA516943; AF073881; LIP-STYX; mMTMH3; MTMR8
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Mtmr9 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 210376). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC046275 , NM_177594 , NM_177594.1 , BC145946 , BC145948
UniProt ID:	Q9Z2D0
Summary:	Acts as an adapter for myotubularin-related phosphatases (PubMed:12890864). Increases lipid phosphatase MTMR6 catalytic activity, specifically towards phosphatidylinositol 3,5-bisphosphate, and MTMR6 binding affinity for phosphorylated phosphatidylinositols (By similarity). Positively regulates lipid phosphatase MTMR7 catalytic activity (PubMed:12890864). The formation of the MTMR6-MTMR9 complex, stabilizes both MTMR6 and MTMR9 protein levels (By similarity). Plays a role in the late stages of macropinocytosis possibly by regulating MTMR6-mediated dephosphorylation of phosphatidylinositol 3-phosphate in membrane ruffles (By similarity). Negatively regulates DNA damage-induced apoptosis, in part via its association with MTMR6 (By similarity). Does not bind mono-, di- and tri-phosphorylated phosphatidylinositols, phosphatidic acid and phosphatidylserine (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).