

Product datasheet for **TL506448**

A3galt2 Mouse shRNA Plasmid (Locus ID 215493)

Product data:

Product Type:	shRNA Plasmids
Product Name:	A3galt2 Mouse shRNA Plasmid (Locus ID 215493)
Locus ID:	215493
Synonyms:	Gm433; iGb3; iGb3S
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	A3galt2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 215493). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001009819 , NM_001009819.1 , NM_001009819.2 , BC140395 , BC146478
UniProt ID:	Q3V1N9
Summary:	Synthesizes the galactose-alpha(1,3)-galactose group on the glycosphingolipid isoglobotrihexosylceramide or isogloboside 3 (iGb3) by catalyzing the transfer of galactose from UDP-Galactose to its acceptor molecule Gal-beta-1,4-Glc-ceramide. Can also catalyze the addition of galactose to iGb3 itself to form polygalactose structures. Synthesis of iGb3 is the initial step in the formation of the isoglobo-series glycolipid pathway and is the precursor to isogloboside 4 (iGb4) and isoForssman glycolipids. Can glycosylate only lipids and not proteins and is solely responsible for initiating the synthesis of isoglobo-series glycosphingolipids.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).