

Product datasheet for TL506569V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Ptpdc1 Mouse shRNA Lentiviral Particle (Locus ID 218232)

Product data:

Product Type: shRNA Lentiviral Particles

Product Name: Ptpdc1 Mouse shRNA Lentiviral Particle (Locus ID 218232)

Locus ID: 218232

Synonyms: Al843923; AW456874; Naa-1; Ptpcd1

Vector: pGFP-C-shLenti (TR30023)

Format: Lentiviral particles

Components: AW456874 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1

scramble control), 0.5 ml each, >10^7 TU/ml.

RefSeq: BC066081, NM 001301781, NM 207232, NM 207232.1, NM 207232.2, NM 001301781.1,

BC068267

UniProt ID: Q6NZK8

Summary: The protein encoded by this gene is a centrosomal mitotic phosphatase. This protein has

been implicated in centrosomal duplication and cytokinesis. In addition, knockdown of expression levels in non-cycling cells results in extra long cilia, suggesting that this protein may function in regulating cilia length independent of a function in cell cycle control. Alternative splicing results in multiple transcript variants encoding different isoforms.

[provided by RefSeq, Sep 2014]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our custom shRNA service.





Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).