

Product datasheet for **TL506599**

Khryn Mouse shRNA Plasmid (Locus ID 219094)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Khryn Mouse shRNA Plasmid (Locus ID 219094)
Locus ID:	219094
Synonyms:	6330555F21; 9130227C08Rik; AI605259; Kiaa0323; mKIAA0323
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Khryn - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 219094). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC116260 , BC116261 , NM_027143 , NR_138134 , NM_027143.1 , NM_027143.2 , NM_027143.3 , BC034330
UniProt ID:	Q80U38
Summary:	This gene encodes a protein with a C-terminal RNA modifying domain that belongs to a family of ribonucleases typified by eukaryotic Nedd4-binding protein1 and the bacterial YacP-like nucleases (NYN). The NYN domain shares a common protein fold with two other groups of nucleases, the PiIT N-terminal nuclease and FLAP nuclease superfamilies. In addition to the NYN domain, the protein encoded by this gene also contains an N-terminal K homology RNA-binding domain. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).