

Product datasheet for **TL506819V**

Asx11 Mouse shRNA Lentiviral Particle (Locus ID 228790)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Asx11 Mouse shRNA Lentiviral Particle (Locus ID 228790)
Locus ID:	228790
Synonyms:	mKIAA0978
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Asx11 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001039939 , NM_001039939.1 , BC026205 , BC059028 , BC076565 , BC172612 , BC172725 , NM_001039939.2
UniProt ID:	P59598
Summary:	Probable Polycomb group (PcG) protein involved in transcriptional regulation mediated by ligand-bound retinoic acid receptors (RARs) and peroxisome proliferator-activated receptor gamma (PPARG). Acts as a coactivator of RARA and RXRA through association with NCOA1. Acts as a corepressor for PPARG and suppresses its adipocyte differentiation-inducing activity. Non-catalytic component of the PR-DUB complex, a complex that specifically mediates deubiquitination of histone H2A monoubiquitinated at 'Lys-119' (H2AK119ub1). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).