

## Product datasheet for **TL506827**

### Arfgap1 Mouse shRNA Plasmid (Locus ID 228998)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Arfgap1 Mouse shRNA Plasmid (Locus ID 228998)
Locus ID:	228998
Synonyms:	AI115377; Arf1gap
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Arfgap1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 228998). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC052922</a> , <a href="#">BC059817</a> , <a href="#">NM_001177706</a> , <a href="#">NM_001177707</a> , <a href="#">NM_001177708</a> , <a href="#">NM_001177709</a> , <a href="#">NM_001177710</a> , <a href="#">NM_145760</a> , <a href="#">NM_145760.1</a> , <a href="#">NM_145760.2</a> , <a href="#">NM_145760.3</a> , <a href="#">NM_001177706.1</a> , <a href="#">NM_001177707.1</a> , <a href="#">NM_001177708.1</a> , <a href="#">NM_001177709.1</a> , <a href="#">NM_001177710.1</a> , <a href="#">BC016228</a> , <a href="#">BC027794</a> , <a href="#">NM_001363035</a> , <a href="#">NM_001363036</a> , <a href="#">NM_001363037</a> , <a href="#">NM_001177706.2</a>
UniProt ID:	<a href="#">Q9EPJ9</a>
Summary:	GTPase-activating protein (GAP) for the ADP ribosylation factor 1 (ARF1). Involved in membrane trafficking and /or vesicle transport. Promotes hydrolysis of the ARF1-bound GTP and thus, is required for the dissociation of coat proteins from Golgi-derived membranes and vesicles, a prerequisite for vesicle's fusion with target compartment. Probably regulates ARF1-mediated transport via its interaction with the KDELR proteins and TMED2. Overexpression induces the redistribution of the entire Golgi complex to the endoplasmic reticulum, as when ARF1 is deactivated. Its activity is stimulated by phosphoinositides and inhibited by phosphatidylcholine (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).