

## Product datasheet for **TL507069V**

### Gpm6a Mouse shRNA Lentiviral Particle (Locus ID 234267)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Gpm6a Mouse shRNA Lentiviral Particle (Locus ID 234267)
Locus ID:	234267
Synonyms:	Gpm6; M6A
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Gpm6a - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">BC023461</a> , <a href="#">BC024759</a> , <a href="#">BC029683</a> , <a href="#">BC029711</a> , <a href="#">BC033357</a> , <a href="#">BC033394</a> , <a href="#">BC043130</a> , <a href="#">NM_001253754</a> , <a href="#">NM_001253756</a> , <a href="#">NM_153581</a> , <a href="#">NM_153581.1</a> , <a href="#">NM_153581.2</a> , <a href="#">NM_153581.3</a> , <a href="#">NM_153581.4</a> , <a href="#">NM_153581.5</a> , <a href="#">NM_153581.6</a> , <a href="#">NM_001253756.1</a> , <a href="#">NM_001253754.1</a>
UniProt ID:	<a href="#">P35802</a>
Summary:	Involved in neuronal differentiation, including differentiation and migration of neuronal stem cells. Plays a role in neuronal plasticity and is involved in neurite and filopodia outgrowth, filopodia motility and probably synapse formation. Gpm6a-induced filopodia formation involves mitogen-activated protein kinase (MAPK) and Src signaling pathways. Conflictingly, PubMed:22162747 reports that induced cellular protrusions are simple membrane-wrapped tubules without actin or tubulin-based cytoskeletons and with Gpm6a gliding along membrane edges indicative for a function in actin-independent membrane deformation. May be involved in neuronal NGF-dependent Ca <sup>2+</sup> influx. May be involved in regulation of endocytosis and intracellular trafficking of G-protein-coupled receptors (GPCRs); enhances internalization and recycling of mu-type opioid receptor.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).