

## Product datasheet for **TL507141**

### Topbp1 Mouse shRNA Plasmid (Locus ID 235559)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Topbp1 Mouse shRNA Plasmid (Locus ID 235559)
Locus ID:	235559
Synonyms:	1110031N14Rik; 2810429C13Rik; AI256758; D430026L04Rik; mKIAA0259
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Topbp1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 235559). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC049797</a> , <a href="#">NM_176979</a> , <a href="#">NM_176979.3</a> , <a href="#">NM_176979.4</a> , <a href="#">NM_176979.5</a> , <a href="#">BC006707</a> , <a href="#">BC007170</a> , <a href="#">BC026608</a> , <a href="#">BC062111</a>
UniProt ID:	<a href="#">Q6ZQF0</a>
Summary:	Required for DNA replication (By similarity). Plays a role in the rescue of stalled replication forks and checkpoint control (PubMed:14718568). Binds double-stranded DNA breaks and nicks as well as single-stranded DNA (By similarity). Recruits the SWI/SNF chromatin remodeling complex to E2F1-responsive promoters. Down-regulates E2F1 activity and inhibits E2F1-dependent apoptosis during G1/S transition and after DNA damage (By similarity). Induces a large increase in the kinase activity of ATR (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).