

Product datasheet for **TL507254**

Abcc4 Mouse shRNA Plasmid (Locus ID 239273)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Abcc4 Mouse shRNA Plasmid (Locus ID 239273)
Locus ID:	239273
Synonyms:	ABCC4-N1; D630049P08Rik; MOATB; MRP4
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Abcc4 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 239273). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001033336 , NM_001163675 , NM_001163676 , NM_001163675.1 , NM_001163676.1 , NM_001033336.1 , NM_001033336.2 , NM_001033336.3 , BC150822 , BC171974
UniProt ID:	Q3TZN9
Summary:	ATP-dependent transporter of the ATP-binding cassette (ABC) family that actively extrudes physiological compounds and xenobiotics from cells. Transports a range of endogenous molecules that have a key role in cellular communication and signaling, including cyclic nucleotides such as cyclic AMP (cAMP) and cyclic GMP (cGMP), bile acids, steroid conjugates, urate, and prostaglandins. Mediates also the ATP-dependent efflux of glutathione conjugates such as leukotriene C4 (LTC4) and leukotriene B4 (LTB4). The presence of GSH is necessary for the ATP-dependent transport of LTB4, whereas GSH is not required for the transport of LTC4. Mediates the cotransport of bile acids with reduced glutathione (GSH). Transports a wide range of drugs and their metabolites, including anticancer, antiviral and antibiotics molecules (Probable). Confers resistance to anticancer agents (Probable).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).