

Product datasheet for **TL507522**

Fktn Mouse shRNA Plasmid (Locus ID 246179)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Fktn Mouse shRNA Plasmid (Locus ID 246179)
Locus ID:	246179
Synonyms:	D830030O17Rik; Fcmd
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Fktn - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 246179). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC017538 , NM_139309 , NM_139309.1 , NM_139309.2 , NM_139309.3 , NM_139309.4 , NM_001363126 , NM_001363127 , NM_001363128
UniProt ID:	Q8R507
Summary:	Catalyzes the transfer of CDP-ribitol to the distal N-acetylgalactosamine of the phosphorylated O-mannosyl trisaccharide (N-acetylgalactosamine-beta-3-N-acetylglucosamine-beta-4-(phosphate-6-)mannose), a carbohydrate structure present in alpha-dystroglycan (DAG1) (PubMed:12471058). This constitutes the first step in the formation of the ribitol 5-phosphate tandem repeat which links the phosphorylated O-mannosyl trisaccharide to the ligand binding moiety composed of repeats of 3-xylosyl-alpha-1,3-glucuronic acid-beta-1 (By similarity). Required for normal location of POMGNT1 in Golgi membranes, and for normal POMGNT1 activity (PubMed:19017726). May interact with and reinforce a large complex encompassing the outside and inside of muscle membranes (PubMed:19017726, PubMed:22922256). Could be involved in brain development (Probable). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).