

Product datasheet for **TL508133V**

Adcy10 Mouse shRNA Lentiviral Particle (Locus ID 271639)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Adcy10 Mouse shRNA Lentiviral Particle (Locus ID 271639)
Locus ID:	271639
Synonyms:	4930431D04Rik; 4931412F17; sAC; Sacy
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Adcy10 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_173029 , NM_001357427 , NM_173029.1 , NM_173029.2 , NM_173029.3 , BC138448 , BC145283
UniProt ID:	Q8C0T9
Summary:	Catalyzes the formation of the signaling molecule cAMP. May function as sensor that mediates responses to changes in cellular bicarbonate and CO ₂ levels (By similarity). Has a critical role in mammalian spermatogenesis by producing the cAMP which regulates cAMP-responsive nuclear factors indispensable for sperm maturation in the epididymis. Induces capacitation, the maturational process that sperm undergo prior to fertilization (PubMed:14976244, PubMed:16054031). Involved in ciliary beat regulation (By similarity). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).