

Product datasheet for **TL508381**

Mical2 Mouse shRNA Plasmid (Locus ID 320878)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Mical2 Mouse shRNA Plasmid (Locus ID 320878)
Locus ID:	320878
Synonyms:	5330438E18Rik; 9530064J02; MICAL-2; mKIAA0750
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Mical2 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 320878). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC111895 , BC112415 , NM_001193305 , NM_177282 , NM_177282.1 , NM_177282.2 , NM_177282.3 , NM_177282.4 , NM_177282.5 , NM_001193305.1 , NM_027587.2
UniProt ID:	Q8BML1
Summary:	Nuclear monooxygenase that promotes depolymerization of F-actin by mediating oxidation of specific methionine residues on actin to form methionine-sulfoxide, resulting in actin filament disassembly and preventing repolymerization (PubMed:23911929, PubMed:23927065). In the absence of actin, it also functions as a NADPH oxidase producing H ₂ O ₂ (By similarity). Acts as a key regulator of the SRF signaling pathway elicited by nerve growth factor and serum: mediates oxidation and subsequent depolymerization of nuclear actin, leading to increase MKL1/MRTF-A presence in the nucleus and promote SRF:MKL1/MRTF-A-dependent gene transcription. Does not activate SRF:MKL1/MRTF-A through RhoA (By similarity). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:

GFP signal was observed under microscope at 48 hours after transduction of [TL508381C] virus into HEK293 cells. [TL508381C] virus was prepared using lenti-shRNA [TL508381C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL508381D] virus into HEK293 cells. [TL508381D] virus was prepared using lenti-shRNA [TL508381D] and [TR30037] packaging kit.