

## Product datasheet for **TL508431**

### Dennd1b Mouse shRNA Plasmid (Locus ID 329260)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Dennd1b Mouse shRNA Plasmid (Locus ID 329260)
Locus ID:	329260
Synonyms:	4632404N19Rik; 4930467M19Rik; 6820401H01Rik; F730008N07Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Dennd1b - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 329260). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC042698</a> , <a href="#">NM_001166501</a> , <a href="#">NM_181347</a> , <a href="#">NM_183217</a> , <a href="#">NM_181347.1</a> , <a href="#">NM_181347.2</a> , <a href="#">NM_181347.3</a> , <a href="#">NM_001166501.1</a> , <a href="#">BC141389</a> , <a href="#">BC141388</a>
UniProt ID:	<a href="#">Q3U1T9</a>
Summary:	Guanine nucleotide exchange factor (GEF) for RAB35 that acts as a regulator of T-cell receptor (TCR) internalization in TH2 cells. Acts by promoting the exchange of GDP to GTP, converting inactive GDP-bound RAB35 into its active GTP-bound form. Plays a role in clathrin-mediated endocytosis. Controls cytokine production in TH2 lymphocytes by controlling the rate of TCR internalization and routing to endosomes: acts by mediating clathrin-mediated endocytosis of TCR via its interaction with the adapter protein complex 2 (AP-2) and GEF activity. Dysregulation leads to impaired TCR down-modulation and recycling, affecting cytokine production in TH2 cells.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).