

## Product datasheet for **TL508522**

### Fa2h Mouse shRNA Plasmid (Locus ID 338521)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Fa2h Mouse shRNA Plasmid (Locus ID 338521)
Locus ID:	338521
Synonyms:	FAAH; Faxdc1; G630055L08Rik
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Fa2h - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 338521). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC128080</a> , <a href="#">BC128081</a> , <a href="#">NM_178086</a> , <a href="#">NM_178086.1</a> , <a href="#">NM_178086.2</a> , <a href="#">NM_178086.3</a> , <a href="#">BC026400</a> , <a href="#">BC026629</a> , <a href="#">BC046985</a> , <a href="#">BC111912</a>
UniProt ID:	<a href="#">Q5MPPQ</a>
Summary:	Catalyzes stereospecific hydroxylation of free fatty acids at the C-2 position to produce (R)-2-hydroxy fatty acids, which are building blocks of sphingolipids and glycosphingolipids common in neural tissue and epidermis (PubMed:15658937, PubMed:16998236). Plays an essential role in the synthesis of galactosphingolipids of the myelin sheath (PubMed:15658937, PubMed:18815260). Responsible for the synthesis of sphingolipids and glycosphingolipids involved in the formation of epidermal lamellar bodies, critical for skin permeability barrier (By similarity). Participates in the synthesis of glycosphingolipids and a fraction of type II wax diesters in sebaceous gland, specifically regulating hair follicle homeostasis (PubMed:21628453). Involved in the synthesis of sphingolipids of plasma membrane rafts, controlling lipid raft mobility and trafficking of raft-associated proteins (PubMed:22517924).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).