

Product datasheet for **TL508673**

Noto Mouse shRNA Plasmid (Locus ID 384452)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Noto Mouse shRNA Plasmid (Locus ID 384452)
Locus ID:	384452
Synonyms:	Flh; MmNot; Not; tc
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Noto - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 384452). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001007472 , NM_001007472.1 , NM_001007472.2 , BC139155 , BC139154
UniProt ID:	Q5TIS6
Summary:	Transcription factor that controls node morphogenesis (PubMed:15231714, PubMed:17884984, PubMed:18061569, PubMed:22357932). Acts downstream of both FOXA2 and Brachyury (T) during notochord development (PubMed:15231714). Is essential for cilia formation in the posterior notochord (PNC) and for left-right patterning; acts upstream of FOXJ1 and RFX3 in this process and is required for the expression of various components important for axonemal assembly and function (PubMed:17884984). Plays a role in regulating axial versus paraxial cell fate (PubMed:18061569). Activates the transcription of ciliary proteins C11orf97 homolog, FAM183B and SPACA9 in the embryonic ventral node (PubMed:27914912).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).