

Product datasheet for **TL509021**

Fam65b Mouse shRNA Plasmid (Locus ID 193385)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Fam65b Mouse shRNA Plasmid (Locus ID 193385)
Locus ID:	193385
Synonyms:	1700108N18Rik; 6330500D04Rik; AI225904; E430013J17Rik; Fam65b; RP23-461O24.1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Fam65b - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 193385). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001080381 , NM_001286100 , NM_001286101 , NM_029679 , NM_178658 , NM_178658.1 , NM_178658.2 , NM_178658.3 , NM_178658.4 , NM_178658.5 , NM_001080381.1 , NM_029679.1 , NM_029679.2 , NM_001286101.1 , NM_001286100.1 , BC059179 , BC060202
UniProt ID:	Q80U16
Summary:	Acts as an inhibitor of the small GTPase RHOA and plays several roles in the regulation of myoblast and hair cell differentiation, lymphocyte T proliferation and neutrophil polarization (PubMed:25588844, PubMed:27269051). Plays a role in fetal mononuclear myoblast differentiation by promoting filopodia and myotube formation (PubMed:17150207). Maintains naive T lymphocytes in a quiescent state and prevents chemokine-induced T lymphocyte responses, such as cell adhesion, polarization and migration (By similarity). Involved also in the regulation of neutrophil polarization, chemotaxis and adhesion (PubMed:25588844). Required for normal development of inner and outer hair cell stereocilia within the cochlea of the inner ear (PubMed:27269051). Plays a role for maintaining the structural organization of the basal domain of stereocilia (PubMed:27269051). Involved in mechanosensory hair cell function (PubMed:27269051). Required for normal hearing (PubMed:27269051).[UniProtKB/Swiss-Prot Function]


[View online »](#)

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).