

## Product datasheet for **TL509032**

### Cdk6 Mouse shRNA Plasmid (Locus ID 12571)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cdk6 Mouse shRNA Plasmid (Locus ID 12571)
Locus ID:	12571
Synonyms:	AI504062; Crk; Crk2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Cdk6 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 12571). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC119360</a> , <a href="#">BC120595</a> , <a href="#">NM_009873</a> , <a href="#">NM_009873.1</a> , <a href="#">NM_009873.2</a> , <a href="#">NM_009873.3</a> , <a href="#">BC027662</a>
UniProt ID:	<a href="#">Q64261</a>
Summary:	This gene encodes a member of the cyclin dependent kinase family of proteins that play important roles in the progression and regulation of the cell cycle. The encoded protein binds to a D-type cyclin to form an active kinase complex to regulate progression through the G1 phase of the cell cycle. Mice lacking the encoded protein exhibit thymic and splenic hypoplasia, and hematopoietic defects such as reduced number of megakaryocytes and erythrocytes. A pseudogene of this gene has been defined on chromosome 4. [provided by RefSeq, Aug 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).