

Product datasheet for **TL509081**

Aff1 Mouse shRNA Plasmid (Locus ID 17355)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Aff1 Mouse shRNA Plasmid (Locus ID 17355)
Locus ID:	17355
Synonyms:	9630032B01Rik; Af; Af4; AW319193; Mllt; Mllt2h; R; Rob
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Aff1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 17355). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC131665 , NM_001080798 , NM_133919 , NM_133919.2 , NM_133919.3 , NM_133919.4 , NM_001080798.1 , NM_001080798.2 , BC018434 , BC085290
Summary:	This gene encodes a member of the AF4/ lymphoid nuclear protein related to the Fragile X E syndrome (FRAXE) family of proteins, which have been implicated in human childhood lymphoblastic leukemia, fragile chromosome X intellectual disability, and ataxia. It is the prevalent mixed-lineage leukemia fusion gene associated with spontaneous acute lymphoblastic leukemia. Members of this family have three conserved domains: an N-terminal homology domain, an AF4/ lymphoid nuclear protein domain, and a C-terminal homology domain. Knockout of the mouse gene by homologous recombination severely affects early events in lymphopoiesis, including precursor proliferation or recruitment, but is dispensable for terminal differentiation. In addition, an autosomal dominant missense mutation results in several phenotypes including ataxia and adult-onset Purkinje cell loss in the cerebellum, indicating a role in Purkinje cell maintenance and function. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).