

## Product datasheet for **TL509085**

### Clec16a Mouse shRNA Plasmid (Locus ID 74374)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Clec16a Mouse shRNA Plasmid (Locus ID 74374)
Locus ID:	74374
Synonyms:	4932416N17Rik; curt; mKIAA0350
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Clec16a - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 74374). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001204229</a> , <a href="#">NM_177562</a> , <a href="#">NM_177562.1</a> , <a href="#">NM_177562.2</a> , <a href="#">NM_177562.3</a> , <a href="#">NM_177562.4</a> , <a href="#">NM_177562.5</a> , <a href="#">NM_001204229.1</a> , <a href="#">BC138907</a> , <a href="#">BC049280</a> , <a href="#">BC138909</a>
UniProt ID:	<a href="#">Q80U30</a>
Summary:	Regulator of mitophagy through the upstream regulation of the RNF41/NRDP1-PRKN pathway. Mitophagy is a selective form of autophagy necessary for mitochondrial quality control. The RNF41/NRDP1-PRKN pathway regulates autophagosome-lysosome fusion during late mitophagy. May protect RNF41/NRDP1 from proteosomal degradation, RNF41/NRDP1 which regulates proteosomal degradation of PRKN. Plays a key role in beta cells functions by regulating mitophagy/autophagy and mitochondrial health.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).