

## Product datasheet for **TL509101**

### Lyar Mouse shRNA Plasmid (Locus ID 17089)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Lyar Mouse shRNA Plasmid (Locus ID 17089)
Locus ID:	17089
Synonyms:	MLZ-264
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Lyar - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 17089). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">BC116923</a> , <a href="#">BC116927</a> , <a href="#">NM_025281</a> , <a href="#">NM_025281.1</a> , <a href="#">NM_025281.2</a> , <a href="#">NM_025281.3</a>
UniProt ID:	<a href="#">Q08288</a>
Summary:	Plays a role in the maintenance of the appropriate processing of 47S/45S pre-rRNA to 32S/30S pre-rRNAs and their subsequent processing to produce 18S and 28S rRNAs (By similarity). Also acts at the level of transcription regulation. Along with PRMT5, binds embryonic globin promoter (By similarity). Represses the expression of embryonic globin Hbb-y gene (PubMed:25092918). In neuroblastoma cells, may also repress the expression of oxidative stress genes, including CHAC1, HMOX1, SLC7A11, ULBP1 and that encoding the small nucleolar RNA SNORD41 (By similarity). Preferentially binds to a DNA motif containing 5'-GGTTAT-3' (By similarity). Stimulates phagocytosis of photoreceptor outer segments by retinal pigment epithelial cells (PubMed:25735755). Prevents NCL self-cleavage, maintaining a normal steady-state level of NCL protein in undifferentiated embryonic stem cells (ESCs), which in turn is essential for ESC self-renewal (PubMed:19489080).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).