

Product datasheet for TL509117

Parp1 Mouse shRNA Plasmid (Locus ID 11545)

Product data:

Product Type: shRNA Plasmids **Product Name:** Parp1 Mouse shRNA Plasmid (Locus ID 11545) Locus ID: 11545 5830444G22Rik; Adprp; Adprt1; Al893648; ARTD1; C80510; PARP; parp-1; PPOL; sPARP-1 Synonyms: pGFP-C-shLenti (TR30023) Vector: E. coli Selection: Chloramphenicol (34 ug/ml) Mammalian Cell Puromycin Selection: Format: Lentiviral plasmids Parp1 - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 11545). **Components:** 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free. BC012041, NM 007415, NM 007415.1, NM 007415.2, NM 007415.3 RefSeq: shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our custom shRNA service. Performance OriGene guarantees that the sequences in the shRNA expression cassettes are verified to Guaranteed: correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over gPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples. For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn