

Product datasheet for **TL509164**

Xiap Mouse shRNA Plasmid (Locus ID 11798)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Xiap Mouse shRNA Plasmid (Locus ID 11798)
Locus ID:	11798
Synonyms:	1110015C02Rik; A; Aipa; Api3; Bir; Birc4; I; IAP3; IL; ILP-1; MIHA
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	Xiap - Mouse, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 11798). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001301639 , NM_001301641 , NM_009688 , NR_125870 , NM_009688.1 , NM_009688.2 , NM_009688.3 , NM_001301639.1 , NM_001301641.1 , BC138528 , BC032190 , BC145861
UniProt ID:	Q60989
Summary:	The protein encoded by this gene is a member of the inhibitor of apoptosis (IAP) family of proteins. While first identified for its role in blocking apoptosis, this protein modulates many other signaling processes including nuclear factor kappa-light-chain-enhancer of activated B cells (NF-κB) pathways and inflammatory responses. This protein blocks apoptosis by binding and inhibiting target caspases after they have been activated. Binding occurs to some, but not all, caspases. This protein has several conserved regions, including baculoviral IAP repeat (BIR) motifs and a RING finger E3 ligase domain. In humans, mutations in this gene are linked to immunodeficiency in X-linked lymphoproliferative syndrome type-2 (XLP-2). A pseudogene of this gene is found on chromosome 7. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).