

Product datasheet for **TL509204V**

Nlgn1 Mouse shRNA Lentiviral Particle (Locus ID 192167)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Nlgn1 Mouse shRNA Lentiviral Particle (Locus ID 192167)
Locus ID:	192167
Synonyms:	6330415N05Rik; BB179718; mKIAA1070; NL1; Nlg1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	Nlgn1 - Mouse shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC098461 , NM_001163387 , NM_138666 , NM_001357095 , NM_001357096 , NM_001357097 , NM_138666.1 , NM_138666.2 , NM_138666.3 , NM_001163387.1 , NM_001163387.2 , BC005523 , NM_138666.4
UniProt ID:	Q99K10
Summary:	Cell surface protein involved in cell-cell-interactions via its interactions with neurexin family members. Plays a role in synapse function and synaptic signal transmission, and probably mediates its effects by recruiting and clustering other synaptic proteins. May promote the initial formation of synapses, but is not essential for this. In vitro, triggers the de novo formation of presynaptic structures. May be involved in specification of excitatory synapses. Required to maintain wakefulness quality and normal synchrony of cerebral cortex activity during wakefulness and sleep (PubMed:23716671).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).